

Vitreia[®] Enterprise Suite 1.3

Vitreia Enterprise Suite 1.3 delivers new clinical capabilities and an enriched set of workflow tools that build on the strength of our reliable, easy-to-use thin client solution. The complete set of clinical solutions is available across your enterprise - anywhere, anytime - to support efficient, effective diagnosis and improve patient outcomes.

Key Features

New CT Cardiac Multi-Chamber Analysis Applications

- Automated atrial and ventricular chamber segmentation and functional measurements
- Leverages low-dose CT scanning protocols

Cardiovascular CTA Analysis

- Vessel Probe & Vessel Walk enhancements ease manipulation and navigation of complicated cardiac and vascular anatomy
- Multi-point vessel seed points, trimming and editing increase read flow accuracy and speed

Workflow Tools to Boost your Productivity

- Single-click Volume Segmentation
- Web-based access to the MeVis Visia[™] CT Lung CAD and Medicsight ColonCAD[™] results within Vitrea Enterprise Suite (ColonCAD[™] available outside the U.S. only)
- Edge detecting contouring and automated table removal speed Multi-Region Segmentation read flow by anticipating your needs

PET/CT Enhancements

- 3D Navigation/triangulation speeds analysis of lesions for SUV characterization

Data Manager – “Save As”

- Save snapshots and movies to your desktop – anywhere

DICOM Print

- Quick and easy creation of DICOM films

Additional Platform Support

- Microsoft[®] Windows[®] 7 clients
- Support for SQL 2008 and SQL Express editions