

IHE Integration Statement			
Vendor	Product Name	Version	Date
Vital Images, Inc	Vitrea	5.1	1/18/2010
Vital Images, Inc	fX	3.0	1/18/2010
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Access to Radiology Information	Image Display	Multiple Sources	
Portable Data for Imaging	Image Display	None	
	Portable Media Creator	None	
	Portable Media Importer	None	
Scheduled Workflow	Evidence Creator	None	
	Image Display	None	
Teaching File and Clinical Trial Export	Export Selector	None	
	Receiver	None	
Internet address for vendor's IHE information: http://www.vitalimages.com/Our_Company/Compliance_Documents.aspx			
Link to Standards Conformance Statements for the Implementation			
DICOM	http://www.vitalimages.com/Our_Company/Compliance_Documents.aspx		
Links to General Information on IHE			
In North America: http://www.ihe.net	In Europe: http://www.ihe-europe.org		In Japan: http://www.jira-net.or.jp/ihe-j

End of Document