



## Vital Life-Cycle Document

Vendor	Product Name	Version	Date
Vital Images, Inc.	Vitrea2	4.1	11/27/2007
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			

Integration Profiles Implemented	Actors Implemented	Options Implemented
Access to Radiology Information	Image Display	Multiple Sources
Portable Data for Imaging	Image Display	None
	Portable Media Creator	None
	Portable Media Importer	None
Scheduled Workflow	Evidence Creator	None
	Image Display	None
Teaching File and Clinical Trial Export	Export Selector	None
	Receiver	None
<b>Internet address for vendor's IHE information:</b> <a href="http://www.vitalimages.com/Solutions/IHE_Integration.aspx">http://www.vitalimages.com/Solutions/IHE_Integration.aspx</a>		
<b>Link to Standards Conformance Statements for the Implementation</b>		

---

End of Document